


## Nurse Call System

Main Receiver / Central Controller /  
Bed Console / Operation Room Interphone /  
Radiation Room Interphone



### MAIN PANEL

JNS-PSM

- ▶ Power supply and signal relay
- ▶ Automatic UPS connection in case of a blackout



### NURSE STATION

DIGITAL

JNS-70MN

Up To  
**96**  
Beds



- ▶ 7" full touch screen with user-friendly GUI.
- ▶ A simple memo function
- ▶ Group announcement function
- ▶ Hands-free function
- ▶ Installation flexibility [Desktop/Surface-mounted type]

DIGITAL

JNS-4CM



- ▶ Easy put the programming of BED and room number
- ▶ Max 48 channel covered
- ▶ Various and clear status indication by 16 digit 2 lines LED display
- ▶ Desk top

ANALOG

JNS-12/24/36/360



- ▶ LED status indication
- ▶ Hands-set
- ▶ Volume control enabled
- ▶ Channel select function
- ▶ Group announcement function
- ▶ Desk top

# BED HEAD CONSOLE UNIT

## JNS-3000



- Calling nurse station
- Various ways of emergency paging and calling
- Interlocking with corridor lights (CL-302C and others)
- Multi-device featured console (AC outlet, light, IV holder, name card)

## JNS-1000C



- Calling nurse station
- Various ways of emergency paging and calling
- Interlocking with corridor lights (CL-302C and others)
- Multi-device featured console (AC outlet, light)

## JNS-101P/ JNS-4CS



- Paging and calling nurse's waiting room
- Emergency paging and calling
- Interlocking with corridor lights (CL-302C and others)

## CC-200



- Emergency call cord

# EMERGENCY UNIT

## CL-302C



- Corridor light
- 2 Colors mode (Red/Green, Red/White)
- LED

## CL-302I



- Corridor light
- 2 Colors mode (Red/Green)
- LED

## ES-400



- Emergency switch (restroom)
- Push Button Type
- Waterproof
- Flush Mount Type

## PB-500



- Presence switch
- Flush Mount Type

## ES-410



- Emergency switch (restroom)
- Waterproof
- Flush Mount Type

## ES-420



- Emergency switch (restroom)
- Waterproof

## PS-100



- Pillow speaker
- Remote control of On/Off (2 lamps)

# OPERATION ROOM UNIT

## OPERATING ROOM INTERPHONE JNS-360

Designed for operating surgeon or nurse to call to waiting room or other operating room for assistance without using hands.  
Two way communication from operating room to another group.

### Master Station

- LED status indication
- Hands-set
- Volume control enabled
- Channel select function
- Group announcement function
- Desktop
- Dimension (mm): 195 (W) x 117 (H) x 16 (D)




## Sub Station I : JNS-101P



## Sub Station II : OPERATING ROOM CALL SYSTEM JNS-6KO

- A simple name tag
- LED status indication
- Flush-mounted type
- Hands-free function
- Dimension (mm) : 305 x 117 x 16



## OPERATING ROOM ARM SWITCH AS-100

- Operation room
- Hands-free intercom
- Arm switch
- Flush-mounted type
- Dimension (mm) : 195 x 117 x 16



## OPERATING ROOM FOOT SWITCH FS-600

- LED status indication
- Foot switch
- Arm switch
- Dimension (mm) : 573 x 117 x 16 (66)



# SYSTEM DIAGRAM



## SOFTWARE NURSE CALL VIEWER PROGRAM

- The web-based viewer program that provides the records of a patient and an emergency call through a PC in real time
- Software: Google, Fire Fox, iOS, IE
- Connectivity : JNS-70MN
- Available Microsoft Office Excel



0857 3325 3456  
0852 2000 9448



www.bnmedika.com



PT. BINTANG NUSANTARA MEDIKA

Pastikan Anda Mendapat Kartu Garansi Dari