















| Product types 产品类别 | Photo 图片 | P/N Example 型号示例 | P/N Guide 型号识别 |
|--|---|---------------------|--|
| SpO2 Adapter Cables 血氧转接线 |  | P 0 2 21B | P: SpO2 0: SpO2 Adapter Cables 2: 2.2m (1: 1.0m; 3: 3.0m) 21B: NIHON KOHDEN |
| Warranty: one-year warranty for reusable sensors and cables Package: 1pcs/bag Length: 2.2m 保修: 一年 包装: 1个/包 标准长度: 2.2m | | | |

SpO2 Adapter Cables/血氧转接线

| Monitor 监护仪厂家 | Medke P/N 曼克型号 | Use with sensor 适用于 | Monitor side conn. P/N 监护仪端接头 | Sensor side conn. P/N 探头端接头 |
|-------------------|-------------------|------------------------|--|---|
| Artema-S&W 阿提玛 | P0201 | Nellcor non-oximax |  CN01 |  CN49 |
| | P0201A | Ohmeda |  CN02 | |
| Baxter 百特 | P0202 | Nellcor non-oximax/BCI |  CN03 |  CN49 |
| Biolight 宝莱特 | P0205 | Nellcor non-oximax |  CP5 |  CN49 |
| | P0205A | Digital |  CP5 |  CN49 |
| | P0205B | Nellcor Oximax |  CP7 |  CN51 |
| | P0205D | Digital |  CN61 |  CN49 |

| Monitor 监护仪厂家 | Medke P/N 曼克型号 | Use with sensor 适用于 | Monitor side conn. P/N 监护仪端接头 | Sensor side conn. P/N 探头端接头 |
|------------------|-------------------|------------------------|---|--|
| Biolight 宝莱特 | P0205E | BLT Db9 digital |  CN62 |  CN49 |
| | P0205H | Digital |  CN82 | |
| | P0205J | BLT new oximax digital |  CN82 | |
| | P0205L | BLT new oximax analog |  CN82B | |
| | P0205F | Nellcor Oximax |  CN61 | |
| Bionet | P0206 | Nellcor non-oximax |  CP7 |  CN49 |
| BCI-Smith | P0203 | BCI |  CP7 |  CN50 |
| | P0203A 3311 | BCI |  CN81 | |
| Biosys | P0204 | Nellcor non-oximax/BCI |  CN81 |  CN49 |
| Comen 科曼 | P0265D | DB9 digit sensor |  CN93 |  CN49 |
| | P0261M-N | Masim MLNCS |  CN93 |  CN58 |
| Contec 康泰 | P02132 | DB9 digit sensor |  CN452 |  CN49 |
| Co-linn 柯林 | P02131 | Nellcor non-oximax |  CN93 |  CN49 |

SpO2 Adapter Cables/血氧转接线

SpO2 Adapter Cables/血氧转接线

SpO2

| Monitor 监护仪厂家 | Medke P/N 曼克型号 | Use with sensor 适用于 | Monitor side conn. P/N 监护仪端接头 | Sensor side conn. P/N 探头端接头 |
|------------------|------------------------|------------------------|--|---|
| Corplus | P0215U | Masimo MLNCS |  CN490B |  CN49 |
| | P0215U-H | Masimo MLNCS | |  CN58 |
| CSI | P0207 518LD | Nellcor non-oximax/CSI |  CM5 |  CN49 |
| | P0207A 518DD | Nellcor non-oximax/CSI |  CN06 |  CN50 |
| Drager 德尔格 | P0209 5720084 | Nellcor non-oximax |  CM8 | |
| | P0209A 4533840 | Nellcor non-oximax |  CN09 |  CN49 |
| | P0209B 3375834 3368433 | Nellcor non-oximax |  CN10 | |
| | P0209DM MS17522 | MASIM LNCS | |  CN49 |
| | P0209F | Masim new tech |  CN87 | |
| | P0209E | Nellcor oximax | |  CN51 |
| | P0209C MS17041 | MASIM LNOP |  CN87 |  CN55 |
| GE-cormetrics 科洛 | P0210 4033CAX | Nellcor non-oximax |  CN11 |  CN49 |
| GE-Trusignal | P0210KT TS-G3 | GE trusignal |  CN15 |  CN57 |

| Monitor 监护仪厂家 | Medke P/N 曼克型号 | Use with sensor 适用于 | Monitor side conn. P/N 监护仪端接头 | Sensor side conn. P/N 探头端接头 |
|--------------------|--------------------------------|---------------------|--|---|
| GE-Trusignal | P0210MT TS-M3 | GE trusignal |  CN18 |  CN57 |
| | P0210QS TS-H3 | GE trusignal |  CN16 | |
| | P0210PS TS-N3 | GE trusignal |  CN12 |  CN56 |
| GE-Datex Ohmeda 德恩 | P0210B OXY-C3 | Datex sensor |  CN12 |  CN49 |
| | P0210P OXY-SL3 | OXY-TIP | |  CN54 |
| GE- Hellige 海力格 | P0210D 30344358 | Nellcor non-oximax |  CN13 |  CN49 |
| GE-Marquette 马奎 | P0210E 2006644-001 E9004GE | Nellcor non-oximax |  CN14 |  CN49 |
| | P0210F 2021406-001 | Nellcor oximax | | |
| | P0210G Masim tech | Nellcor non-oximax | | |
| | P0210GM 2027263-002 (REF 2016) | Masim LNCS |  CN15 |  CN49 |
| | P0210W | Masimo DB9 | | |
| | P0210WM | Masim New LNCS | | |
| | P0210W-S | Masimo RD SET |  CN15 |  CN448 |

SpO2 Adapter Cables/血氧转接线

| Monitor 监护仪厂家 | Medke P/N 曼克型号 | Use with sensor 适用于 | Monitor side conn. P/N 监护仪端接头 | Sensor side conn. P/N 探头端接头 |
|--|-----------------------|------------------------|---|--|
| GE-Marquette 马奎 | P0210S PC-12-GE | Masim LNOP |  CN15 |  CN55 |
| GE-Ohmeda 欧美达 | P0210H | Ohmeda |  CN16 |  CN49 |
| | P0210K OXY-ES3 | OXY-TIP |  CN15 |  CN54 |
| | P0210M OXY-MC3 | OXY-TIP |  CN18 | |
| | P0210Q OXY-OL3 | OXY-TIP |  CN16 | |
| Goldway, Pace Tech, Generra 金科威 | P0211 | Nellcor non-oximax/BCI |  CM5 |  CN49 |
| | P0211A | Nellcor non-oximax/BCI |  CP5 | |
| | P0211B | Nellcor non-oximax/BCI |  CP7 | |
| | P0211C | Nellcor non-oximax/BCI |  CN19 | |
| Invivo 英威 | P0212 | Nellcor non-oximax |  CM7 |  CN49 |
| | P0212A Masim tech | Nellcor non-oximax | | |
| | P0212AM Masim tech | Masim LNCS | | |
| Kontron 康强 | P0213 | Nellcor non-oximax/BCI |  CN20 |  CN49 |

SpO2 Adapter Cables/血氧转接线



| Monitor 监护仪厂家 | Medke P/N 曼克型号 | Use with sensor 适用于 | Monitor side conn. P/N 监护仪端接头 | Sensor side conn. P/N 探头端接头 |
|----------------------|--------------------------|---|---|---|
| Kontron 康强 | P0213A | Nellcor non-oximax |  CP7T |  CN49 |
| Mindray 迈瑞 | P0218 0515-30-11221 | Nellcor non-oximax/BCI |  CN39 |  CN49 |
| | P0218A Masim tech | Nellcor non-oximax/BCI |  CN39 | |
| | P0218AM 0010-30-12452 | Masim LNCS |  CM5 | |
| | P0218B 512A-30-06074 | Nellcor non-oximax/BCI | | |
| | P0218C | Nellcor non-oximax/BCI |  CP5 | |
| | P0218D | Nellcor non-oximax/BCI |  CP6 | |
| | P0218E 0010-20-42594 | Nellcor non-oximax/BCI |  CP6R | |
| | P0218F Masim tech | Nellcor non-oximax/BCI |  CP6R | |
| | P0218FM 0010-30-42602 | Masim LNCS |  CP6S | |
| | P0218GM 0010-30-42625 | Masim LNCS | | |
| | P0218H 0010-20-42710 | Nellcor non-oximax/BCI |  CN27 | |
| P0218J Masim tech | Nellcor non-oximax/BCI |  CN28 | | |

SpO2

SpO2 Adapter Cables/血氧转接线

SpO2 Adapter Cables/血氧转接线

SpO2

| Monitor 监护仪厂家 | Medke P/N 型号 | Use with sensor 适用于 | Monitor side conn. P/N 监护仪端接头 | Sensor side conn. P/N 探头端接头 |
|---------------|-----------------------------------|------------------------|---|---|
| Mindray 迈瑞 | P0218JM 0010-30-42738 | Masim LNCS |  CN28 | |
| | P0218P 115-020768-00 | Masim LNCS |  CN91 |  CN49 |
| | P0218L | Nellcor oximax |  CP6T | |
| | P0218K 0010-20-42712 Oximax | Nellcor oximax |  CN29 |  CN51 |
| | P0218M | Nellcor oximax |  CP6R |  CN51 |
| | P0218N 9200-30-10707 | Masim Lnop |  CP6R |  CN55 |
| | P0218P-S | Masim RD |  CN91 |  CN448 |
| Masim 玛思莫 | P0215 Masim tech | Nellcor non-oximax/BCI | |  CN49 |
| | P0215B Masim 1005/PC08 | Masim LNOP |  CN22 |  CN55 |
| | P0215L Masim MLNCS | Masim new LNCS | |  CN49 |
| | P0215N 2525 (M-LNC-10) | Masim MLNCS | |  CN58 |
| | P0215G 2059/RED LNC-10 | Masim LNCS |  CN80 |  CN49 |
| | P0215H | Masim MLNCS | |  CN58 |

| Monitor 监护仪厂家 | Medke P/N 型号 | Use with sensor 适用于 | Monitor side conn. P/N 监护仪端接头 | Sensor side conn. P/N 探头端接头 |
|---------------|---------------------------------|---------------------|---|---|
| Masim 玛思莫 | P0215G-S RD rainbow SET Md20 | Masim RD SET |  CN80 | |
| | P0227E-S | Masim RD |  CN69 |  CN448 |
| | P0215AS-S RAD-G | Masim RD SET |  CN443A | |
| Masim 玛思莫 | P0215Q RED 25 LNC | Masim LNCS |  CN441 |  CN49 |
| | P0215Q-S | Masim RD SET | |  CN448 |
| | P0215Q-B RED 25 LNC | Masim LNCS |  CN441B |  CN49 |
| | P0215Q-BS | Masim RD SET | |  CN448 |
| | P0215T | Masim RD Set Rad-97 | |  CN448 |
| | P0215T-B | Masim LNCS |  CN457 |  CN49 |
| Mennen | P0215T-H | Masim MLNCS | |  CN58 |
| | P0217 | Nellcor non-oximax |  CM10 | |
| | P0217A | Nellcor non-oximax |  CN25 |  CN49 |
| | P0217B | Nellcor non-oximax |  CN26 | |

SpO2 Adapter Cables/血氧转接线

| Monitor 监护仪厂家 | Medka P/N 麦克型号 | Use with sensor 适用于 | Monitor side conn. P/N 监护仪端接头 | Sensor side conn. P/N 探头端接头 |
|-------------------|------------------------------------|------------------------|---|---|
| Mennen | P02178M | Masim LNCS |  |  |
| | P0217C | Masim LNOP | CN26 |  |
| Mindray-Datascope | P0208 0012-00-0516-02 | Datascope |  |  |
| | P0215B 0012-00-1099-01 | Masim Lnop |  |  |
| MEK | P0216A | DB9 MEK sensor |  |  |
| | P0216B | Nellcor non-oximax/BCI | CN65 | |
| | P0216C | MEK |  | CN66 |
| Nellcor | P0219 DEC-8/EC-8/ DEC-4/EC-4 | Nellcor non-oximax |  |  |
| | P0019A | Nellcor non-oximax/BCI |  |  |
| | P0219A Oximax cable | Nellcor non-oximax/BCI |  |  |
| | P0219B SCP-10/MC-10 | Nellcor non-oximax |  |  |
| | P0319C DOC-10 | Nellcor oximax |  |  |
| Nihon Kohden 光电 | P0221 | DB9 N.K sensor |  |  |

SpO2 Adapter Cables/血氧转接线

| Monitor 监护仪厂家 | Medka P/N 麦克型号 | Use with sensor 适用于 | Connector 接头型号 | Monitor side conn. P/N 监护仪端接头 | Sensor side conn. P/N 探头端接头 |
|-----------------|-----------------------------------|-----------------------|---|---|---|
| Nihon Kohden 光电 | P0221B JL-900P | DB9 N.K sensor |  |  |  |
| | P0221C | DB9 N.K sensor |  |  | |
| | P0221F | DB9 N.K sensor |  | CP9R | |
| | P0221E | Masim |  |  | |
| | P0221D | Nellcor Oximax | |  | |
| | P0221E-S | Masim RD |  | CN448 | |
| Nonin 燕牌 | P0222 | Nonin sensor | |  | CM6 |
| | P0222A | Nonin sensor |  | CN35 | |
| | P0222B | Nonin sensor |  | CN83, Direct (直连), 3m | |
| Novamatrix 诺匿 | P0223 | Novamatrix sensor |  |  |  |
| Philips, HP 飞利浦 | P0225 M1900B | Nellcor non-oximax |  |  |  |
| | P0225A M1941A | Philips 8 pins sensor |  |  |  |
| | P0225F 1847/LNOP MPI2 Masim | Masim LNOP |  |  |  |

SpO2

SpO2 Adapter Cables/血氧转接线

| Monitor 监护仪厂家 | Medke P/N 曼克型号 | Use with sensor 适用于 | Monitor side conn. P/N 监护仪端接头 | Sensor side conn. P/N 探头端接头 |
|--------------------------------|-----------------------------------|------------------------|----------------------------------|--------------------------------|
| Philips, HP 飞利浦 | P0225B | Nellcor non-oximax | | |
| | P0225BM 2281/LNC MP10 Masim | Masim LNCS sensor | | |
| | P0225L-H | Masim MLNCS | | |
| | P0225L-S | Masim RD Set | | |
| | P0225C M1940A | Philips 8 pins sensor | | |
| | P0225D M1943A/ M1943AL | Nellcor non-oximax | | |
| | P0225E M1943NL | Nellcor oximax | | |
| | P0225E-B | Nellcor oximax | | |
| Solaris-New Tech 索莱瑞-纽康克 | P0220 | DB9 digital sensor | | |
| | P0220A | DB9 digital sensor | | |
| | P0220D | DB9 digital sensor | | |
| Schiller 席勒 | P0226 | BCI sensor | | |
| | P0226A | Nellcor Non-Oximax | | |

SpO2 Adapter Cables/血氧转接线

| Monitor 监护仪厂家 | Medke P/N 曼克型号 | Use with sensor 适用于 | Monitor side conn. P/N 监护仪端接头 | Sensor side conn. P/N 探头端接头 |
|------------------|------------------------|------------------------|----------------------------------|--------------------------------|
| Schiller 席勒 | P0226AM | Masim LNCS sensor | | |
| | P0226B | BCI sensor | | |
| Spacelabs 太空 | P0227 700-0002-00 | Nellcor non-oximax | | |
| | P0227A 700-0030-00 | Nellcor non-oximax | | |
| | P0227B | Nellcor non-oximax | | |
| | P0227BM 700-0906-00 | Masim LNCS | | |
| | P0227D | Nellcor oximax | | |
| | P0227C | Masim LNOP | | |

SpO2