




Reusable SpO2 Sensors  
重复性血氧探头

- Finger clip, Soft tip, Wrap, Multi-Y, Ear clip types available  
指夹, 指套, 包裹, 多功能 Y 型, 耳夹可用
- Latex-free  
不含乳胶
- Available in adult, pediatric, infant and neonate sizes  
成人, 儿童, 婴儿, 新生儿可用
- Different types are Compatible with most of patient monitors  
多款样式的探头兼容大多数病人监护仪
- High performance/price ratio for Reusable SpO2 sensor to save cost  
高性价比的重复性血氧探头, 节约成本
- Easy to clean and repair  
容易清洗和维修



SpO2

Reusable SpO2 Sensors  
重复性血氧探头

Certification  
产品认证


 CE Mark  
 No.: G1 13 11 85432 002  
 FDA 510(k) clearance No.:K152390  
 NMPA 医疗器械注册证 编号:粤械注准20152071060  
 Free Sales Certificate(FSC) No.:20201424  
 自由销售证书 编号:粤食药监械出20201424号

Compliance with  
符合标准

ISO 10993-1:2018  
 ISO 10993-5:2009  
 ISO 10993-10:2010  
 IEC 62366-1:2015  
 EN 60601-1-2:2006/A1:2013  
 EN ISO 80601-2-56:2017/A1:2020

Product types 产品类别	Photo 图片	P/N Example 型号示例	P/N Guide 型号识别
Reusable SpO2 Sensors 重复性血氧探头		P 8 3 25A	P: SpO2 8: Adult soft (9: Adult clip; 7: Pediatric clip; 6: Pediatric soft; 5: Neonate wrap; 4: Multi Y; 3: Ear clip) 3: 3.0m (1: 1.0m) 25A: Philips 8 pins (25: Philips 12 pins; 19: Nellcor 7 pins;19A: Nellcor 9 pins)
	 Adult Ear clip 成人耳夹  Multi Y 多功能Y型  Neonate wrap 新生儿包裹  Pediatric soft 儿童指套  Pediatric clip 儿童指夹  Adult soft 成人指套  Adult clip 成人指夹		
Warranty: one-year warranty for reusable sensors and cables Package: 1pcs/bag 保修: 一年 包装: 1个/包			

Reusable SpO2 Sensors/重复性血氧探头

Monitor 监护仪厂家	Adult dip 成人指夹	Adult soft 成人软指夹	Pedia. Clip 儿童指夹	Pedia. Soft 儿童软指夹	Neonate wrap 新生儿包裹	Multi-Y 多功能Y型	Adult Ear clip 成人耳夹	Connector/Adapter/Len. 接头/转接线/长度
<b>Artema 阿捷玛</b>								
Diascope	P9301	P8301	P7301	P6301	P5301	P4301	P3301	 CN01, Direct (直连), 3m
S&W	P9301A S1033057	P8301A	P7301A	P6301A	P5301A	P4301A	P3301A	 CN02, Direct (直连), 3m
<b>Artemstrong Medica</b>								
AD-900, 1000, 2000, 3000, 4000(BCI tech)	P9303A	P8303A	P7303A	P6303A	P5303A	P4303A	P3303A	 CN81, Direct (直连), 3m
<b>Baxter/Simed 百特</b>								
S-50, S100, S-100e	P9302 600015	P8302	P7302	P6302	P5302	P4302	P3302	 CN03, Direct (直连), 3m
Monarch/Explorer	P9319	P8319	P7319	P6319	P5319	P4319	P3319	 CN30, Direct (直连), 3m
<b>BCI/Smiths</b>								
3100, 3101, 6100, 9100, Advisor, Autocorr, Mini Torr, Mini Torr Plus	P9303	P8303	P7303	P6303	P5303	P4303	P3303	 CP7, Direct (直连), 3m
3300, 3301, 3302, 3303, 3304, 3401, 3180 ChoiceMed MD2000A	P9103A	P8103A	P7103A	P6103A	P5103A	P4103A	P3103A	 CN81, Direct (直连), 1m
	P9303A 3044/3444	P8303A	P7303A	P6303A	P5303A 3025/3026	P4303A	P3303A 3043	 CN81, Direct (直连), 3m
<b>Betterlife 苏州贝莱弗</b>								
Betterlife	P93142	P83142	P73142	P63142	P53142	P43142	P33142	 CN213, Direct (直连), 3m
<b>Bird</b>								
4400	P9307	P8307	P7307	P6307	P5307	P4307	P3307	 CMS, Direct (直连), 3m
<b>Biocare 邦健</b>								
Biocare iM12 (BCI)	P9393	P8393	P7393	P6393	P5393	P4393	P3393	 CN85, Direct (直连), 3m
<b>Biosys</b>								
BPM-200, Vital Sign Monitor SENTRY	P9304	P8304	P7304	P6304	P5304	P4304	P3304	 CN05, Direct (直连), 3m

Reusable SpO2 Sensors/重复性血氧探头

SpO2

Monitor 监护仪厂家	Adult dip 成人指夹	Adult soft 成人软指夹	Pedia. Clip 儿童指夹	Pedia. Soft 儿童软指夹	Neonate wrap 新生儿包裹	Multi-Y 多功能Y型	Adult Ear clip 成人耳夹	Connector/Adapter/Len. 接头/转接线/长度
<b>Bionet</b>								
BM3(before 2005), BM3 Vet	P9306	P8306	P7306	P6306	P5306	P4306	P3306	 CP7-1, Direct (直连), 3m
Bm3	P9306A	P8306A	P7306A	P6306A	P5306A	P4306A	P3306A	 CN63, Direct (直连), 3m
BMS	P9306B	P8306B	P7306B	P6306B	P5306B	P4306B	P3306B	 CN78, Direct (直连), 3m
<b>Biolight 宝莱特</b>								
Biolight M12; LGMD (BCI)	P9305	P8305	P7305	P6305	P5305	P4305	P3305	 CP5, Direct (直连), 3m
Biolight M69; M9000 Biocare Neusoft NSC-M10; Contec CMS9000, CMS6000 (Digital)	P9305A	P8305A	P7305A	P6305A	P5305A	P4305A	P3305A	 CN30, PO205, 1m
M700, M800, Creative SP20 (Digital)	P9105C	P8105C	P7105C	P6105C	P5105C	P4105C	P3105C	 CN30, PO205, 1m
M7000 ,smk-1 , Advanced AM1500 (Digital)	P9305D	P8305D	P7305D	P6305D	P5305D	P4305D	P3305D	
Schiller,Argus,LSM (Nellcor Oximax)	P9305F	P8305F	P7305F	P6305F	P5305F	P4305F	P3305F	 CN61, Direct (直连), 3m
Analogic medical AM 1500 (Digital new) (N)	P9305N	P8305N	P7305N	P6305N	P5305N	P4305N	P3305N	
(Digital)	P9305E	P8305E	P7305E	P6305E	P5305E	P4305E	P3305E	 CN62, Direct (直连), 3m
A8, M9000, V6, Advanced PM-200M, Any View A5 (Digital)	P9305H	P8305H	P7305H	P6305H	P5305H	P4305H	P3305H	
A8 , M9000, V6, Advanced PM-200M, Any View A5, Innomed T12 Plus (Digital new) (N)	P9305J	P8305J	P7305J	P6305J	P5305J	P4305J	P3305J	 CN62, Direct (直连), 3m
A8 (Nellcor Oximax)	P9305G	P8305G	P7305G	P6305G	P5305G	P4305G	P3305G	
AM1500 (Digital new) (N)	P9105K	P8105K	P7105K	P6105K	P5105K	P4105K	P3105K	
Anyview (Analog new) (N)	P9105M	P8105M	P7105M	P6105M	P5105M	P4105M	P3105M	 CN80-C, Direct (直连), 1m

Reusable SpO2 Sensors/重复性血氧探头

Monitor 监护仪厂家	Adult clip 成人指夹	Adult soft 成人软指套	Pedia. Clip 儿童指夹	Pedia. Soft 儿童软指套	Neonate wrap 新生儿包裹	Multi-Y 多功能Y型	Adult ear clip 成人耳夹	Connector/Adapter/Len. 接头/转接线/长度
<b>Biolight 宜莱特</b>								
Anyview S12 (Analog new)	P9305L	P8305L	P7305L	P6305L	P5305L	P4305L	P3305L	 CN2B, Direct (直连), 3m
<b>CAS Medical</b>								
OscilloMate 9002, 9303, 9302S (Nonin tech)	P9322B	P8322B	P7322B	P6322B	P5322B	P4322B	P3322B	 CN83, Direct (直连), 3m
740-750 (Masim tech)	P9315	P8315	P7315	P6315	P5315	P4315	P3315	 CN22, Direct (直连), 3m
<b>Colin 科林</b>								
BP88/306/508/ 608/BX1One (Nellcor tech)	P9319	P8319	P7319	P6319	P5319	P4319	P3319	 CN30, Direct (直连), 3m
BP 88S	P9315	P8315	P7315	P6315	P5315	P4315	P3315	 CN22, Direct (直连), 3m
Bx10ma (Nellcor tech)	P9319B	P8319B	P7319B	P6319B	P5319B	P4319B	P3319B	
<b>Comen 科曼</b>								
Star8000C/8000D, Star 8000, C50 (Digital)	P9361A	P8361A	P7361A	P6361A	P5361A	P4361A	P3361A	 CP6, Direct (直连), 3m
CM8000A, Star-1000(BCI)	P9361	P8361	P7361	P6361	P5361	P4361	P3361	
DOTOP, EDANINS EM9000E, SPRINGS SPR9000G, Jumbor, Hwatime H8, AVEUS AV-Pro(BCI)	P9361B	P8361B	P7361B	P6361B	P5361B	P4361B	P3361B	 CP6R, Direct (直连), 3m
(Digital)	P9361C	P8361C	P7361C	P6361C	P5361C	P4361C	P3361C	
C60 (Nellcor Oximax)	P9361D	P8361D	P7361D	P6361D	P5361D	P4361D	P3361D	 CP6R, Direct (直连), 3m
C60 (Masim tech)	P9361E	P8361E	P7361E	P6361E	P5361E	P4361E	P3361E	
<b>CSI</b>								
501, 503, 504, 507, 508, 581, 1100, 8100, Poet(602-3A), Scholar	P9307 934-10LN	P8307 975AD- 10L	P7307 1563-10L	P6307 975PD- 10L	P5307	P4307 940-10L	P3307 1567-10L	 CM5, Direct (直连), 3m
90 degree angle	P9107A 934SDN	P8107A 975AD-3	P7107A 1563-SD	P6107A 975PD-3	P5107A	P4107A 940SD	P3107A 1567SD	 CN05, See Cable, 1m

Reusable SpO2 Sensors/重复性血氧探头

SpO2

Monitor 监护仪厂家	Adult clip 成人指夹	Adult soft 成人软指套	Pedia. Clip 儿童指夹	Pedia. Soft 儿童软指套	Neonate wrap 新生儿包裹	Multi-Y 多功能Y型	Adult ear clip 成人耳夹	Connector/Adapter/Len. 接头/转接线/长度
<b>Comen 科曼</b>								
STAR 8000A (BCI tech)	P9361J	P8361J	P7361J	P6361J	P5361J	P4361J	P3361J	 CP6S, Direct (直连), 3m
C60 (Masim tech)	P9361K	P8361K	P7361K	P6361K	P5361K	P4361K	P3361K	 CN22, Direct (直连), 3m
C60 (Masim tech)	P9361L	P8361L	P7361L	P6361L	P5361L	P4361L	P3361L	 CM8-1B, Direct (直连), 3m
C50/C60/C80(Masim tech with M* chips)	P9361G	P8361G	P7361G	P6361G	P5361G	P4361G	P3361G	 CN9B, Direct (直连), 3m
Masim	P9361P	P8361P	P7361P	P6361P	P5361P	P4361P	P3361P	
	P9361M	P8361M	P7361M	P6361M	P5361M	P4361M	P3361M	
C100A Telemetry Monitoring (Digital)	P9365D	P8365D	P7365D	P6365D	P5365D	P4365D	P3365D	 CN96, Direct (直连), 3m
C30/C80(Nellcor tech), Star8000 NCS(Nellcor tech)	P9361H	P8361H	P7361H	P6361H	P5361H	P4361H	P3361H	
<b>Contec 康泰</b>								
CMS7000/ 8000/9000 Sinohero H8	P93132	P83132	P73132	P63132	P53132	P43132	P33132	 CN452, Direct (直连), 3m
<b>Choice Electronic 德恩</b>								
All (BCI tech)	P9303A	P8303A	P7303A	P6303A	P5303A	P4303A	P3303A	 CN81, Direct (直连), 3m
M3,MPC300, MMED6000DP	P9342	P8342	P7342	P6342	P5342	P4342	P3342	 CP6, Direct (直连), 3m
Ecomec EPM-50, EPM-80 iplus (Digital)	P9342A	P8342A	P7342A	P6342A	P5342A	P4342A	P3342A	
MD2000B	P9342C	P8342C	P7342C	P6342C	P5342C	P4342C	P3342C	 CN99, Direct (直连), 3m
<b>CO+LINN 柯林</b>								
CO+LINN	P93131	P83131	P73131	P63131	P53131	P43131	P33131	 CN9B, Direct (直连), 3m

Reusable SpO2 Sensors/重复性血氧探头








Monitor 监护仪厂家	Adult clip 成人夹夹	Adult soft 成人软指夹	Pedia. Clip 儿童指夹	Pedia. Soft 儿童软指夹	Neonate wrap 新生儿包裹	Multi-Y 多功能Y型	Adult Ear clip 成人耳夹	Connector/Adapter/Len. 接头/转接线/长度
<b>Charmcare</b>								
Charmcare Cx100 (Nelcor tech)	P9319	P8319	P7319	P6319	P5319	P4319	P3319	 CN30, Direct (直连), 3m
<b>3F (Shanghai)</b>								
Cardioline Vsign 200B(BCI) Biocare IM8	P9354	P8354	P7354	P6354	P5354	P4354	P3354	 CP6, Direct (直连), 3m
<b>Creative 科瑞康</b>								
Creative, Heal force: advance 120, Deluxe 70	P9395	P8395	P7395	P6395	P5395	P4395	P3395	 CP8, Direct (直连), 3m
K12	P93136	P83136	P73136	P63136	P53136	P43136	P33136	 CN475, Direct (直连), 3m
<b>CSI</b>								
503 Spot, 504Plus, 504DX, 506DXN 90 degree angle	P9307A 934-100N	P8307A 975AD-10	P7307A 1563-10D	P6307A 975PD-10	P5307A	P4307A 940-10D	P3307A 1567-10D	 CN06, Direct (直连), 3m
<b>CU Medical</b>								
LiFEGAIN CU-HD1	P93103	P83103	P73103	P63103	P53103	P43103	P33103	 CP8, Direct (直连), 3m
<b>Dolphin 海豚</b>								
2100, 2150 (Digital)	P9162A	P8162A	P7162A	P6162A	P5162A	P4162A	P3162A	 CN7L, See Cable, 1m
	P9362	P8362	P7362	P6362	P5362	P4362	P3362	 CN22, Direct (直连), 3m
<b>Dräger-Siemens 德尔格-西门子</b>								
Dräger Oxipac 2000, 2500 (BCI tech)	P9303A	P8303A	P7303A	P6303A	P5303A	P4303A	P3303A	 CN81, Direct (直连), 3m
Vitalert 1000, Cicero, GammaXL, DeltaXL, Kappa	P9309	P8309	P7309	P6309	P5309	P4309	P3309	 CM8, Direct (直连), 3m
Sirecust, Sirem Spo2	P9309A	P8309A	P7309A	P6309A	P5309A	P4309A	P3309A	 CN09, Direct (直连), 3m
Sc6000, 7000, 9000	P9309B	P8309B	P7309B	P6309B	P5309B	P4309B	P3309B	 CN10, Direct (直连), 3m
(Masim tech)	P9309C	P8309C	P7309C	P6309C	P5309C	P4309C	P3309C	 CN87, Direct (直连), 3m

Reusable SpO2 Sensors/重复性血氧探头

SpO2










Monitor 监护仪厂家	Adult clip 成人夹夹	Adult soft 成人软指夹	Pedia. Clip 儿童指夹	Pedia. Soft 儿童软指夹	Neonate wrap 新生儿包裹	Multi-Y 多功能Y型	Adult Ear clip 成人耳夹	Connector/Adapter/Len. 接头/转接线/长度
<b>Dräger-Siemens 德尔格-西门子</b>								
MS18683 (Nelcor Oximax)	P9309E	P8309E	P7309E	P6309E	P5309E	P4309E	P3309E	 CN87, Direct (直连), 3m
(Masim Oximax)	P9309F	P8309F	P7309F	P6309F	P5309F	P4309F	P3309F	
<b>DST(USA)</b>								
DST(USA)	P9381	P8381	P7381	P6381	P5381	P4381	P3381	 CP6R, Direct (直连), 3m
<b>Edan 理邦</b>								
UTAS, Utas UM300	P9346	P8346	P7346	P6346	P5346	P4346	P3346	 CP8, Direct (直连), 3m
<b>GE-Corometrics 科洛</b>								
Corometrics 188, 129, 129F, 128, 128F	P9319	P8319	P7319	P6319	P5319	P4319	P3319	 CN30, Direct (直连), 3m
Corometrics 120, 510, 511, 118, 556	P9310	P8310	P7310	P6310	P5310	P4310	P3310	 CN11, Direct (直连), 3m
<b>GE-Critikon 强生</b>								
Dinamap 8700/8710/ 8720, Plus, Oxyshuttle	P9310A	P8310A	P7310A	P6310A	P5310A	P4310A	P3310A	 CM7, Direct (直连), 3m
Dinamap 9700 / 9710 / 9720, Plus, Compact, MPS, Select	P9319	P8319	P7319	P6319	P5319	P4319	P3319	 CN30, Direct (直连), 3m
Dinamap ProSeries V1 and Pro 1000 V2 (Nelcor tech)	P9319B	P8319B	P7319B	P6319B	P5319B	P4319B	P3319B	 CN22, Direct (直连), 3m
ProCare; ProSeries V2; Pro 1000 V3 (Masim)	P9315	P8315	P7315	P6315	P5315	P4315	P3315	
Same as above (Nelcor Oximax)	P9119A	P8119A	P7119A	P6119A	P5119A	P4119A	P3119A	 CN31, P0319C, 1m
<b>GE-Datex Ohmeda 德恩</b>								
S/5, AS/3, CS/3 Modular Monitors; with M-series modules; Cardiicap/S, Cardiicap S/5 Light Monitor; Capnomac Ultima Satellite Plus II	P9310B OKY-F4-N	P8310B	P7310B	P6310B	P5310B	P4310B	P3310B OXY-E4-N	 CN12, Direct (直连), 3m
	P9110C OKY-F-DB	P8110C	P7110C	P6110C	P5110C OKY-W-DB	P4110C	P3110C OXY-E-DB	
								 CN04, P0210B, 1m

Reusable SpO2 Sensors/重复性血氧探头

Monitor 监护仪厂家	Adult clip 成人指夹	Adult soft 成人软指套	Pedia. Clip 儿童指夹	Pedia. Soft 儿童软指套	Neonate wrap 新生儿包裹	Multi-Y 多功能Y型	Adult Ear clip 成人耳夹	Connector/Adapter/Len. 接头/转接线/长度
<b>GE-Datex Ohmeda 德恩</b>								
M-NSAT Modules for AS/3, CS/3,S/5 Modular	P9319B	P8319B	P7319B	P6319B	P5319B	P4319B	P3319B	 CN22, Direct (直连), 3m
E-NSAT Modules for i4 Modular Monitors	P9119A	P8119A	P7119A	P6119A	P5119A	P4119A	P3119A	 CN31, P0319C, 3m
<b>GE-Hellige 海力格</b>								
Eagle 1000, 3000, SMU, SMK, SMV, cardioserv	P9310D	P8310D	P7310D	P6310D	P5310D	P4310D	P3310D	 CN13, Direct (直连), 3m
<b>GE-Marquette 马奎</b>								
Dash 2000/3000/4000; Solar8000/8000M, i9500 with Tram Module451/451N/851/851N;CardioLab/Mac-Labwith Tram Module 451/451N;Corometrics 120/2120 series; Eagle 1000/3000/4000;Responder 3000 (Nellcor tech) Oxismart	P9310E	P8310E	P7310E	P6310E	P5310E	P4310E	P3310E	 CN14, Direct (直连), 3m
Dash 3000/4000/5000 (OxiMax);Solar 8000/8000M, 19500 with Tram module 451N5/851N5	P9310F	P8310F	P7310F	P6310F	P5310F	P4310F	P3310F	 CN15, Direct (直连), 3m
Dash 3000/4000/5000 Solar 8000/8000M, i9500 with Tram Module 451M/851M CardioLab /Mac-Lab with Tram Module 451M; Corometrics 259(Masim)	P9310G	P8310G	P7310G	P6310G	P5310G	P4310G	P3310G	
GE CARESCAPE V100 , GE CARESCAPE B650 Masim Tech (with M* chip)	P9310W	P8310W	P7310W	P6310W	P5310W	P4310W	P3310W	
<b>GE-Ohmeda 欧姆达</b>								
Ohmeda TuffSat;3775, 37xx Series, 3800, 3900,3900P,4700 and 5250 respiratory gas monitors;Cardiocap/5	P9310H	P8310H	P7310H	P6310H	P5310H	P4310H	P3310H	 CN16, Direct (直连), 3m
Ohmeda TuffSat	P9110H	P8110H	P7110H	P6110H	P5110H	P4110H	P3110H	 CN16, Direct (直连), 1m

Reusable SpO2 Sensors/重复性血氧探头

SpO2






Monitor 监护仪厂家	Adult clip 成人指夹	Adult soft 成人软指套	Pedia. Clip 儿童指夹	Pedia. Soft 儿童软指套	Neonate wrap 新生儿包裹	Multi-Y 多功能Y型	Adult Ear clip 成人耳夹	Connector/Adapter/Len. 接头/转接线/长度
<b>GE-Ohmeda 欧姆达</b>								
Ohmeda TuffSat;3775, 37xx Series, 3800, 3900,3900P,4700 and 5250 respiratory gas monitors;Cardiocap/5	P9110J	P8110J	P7110J	P6110J	P5110J	P4110J	P3110J	 CN04, P0210H, 1m
Ohmeda TruSat, Ohmeda TruSat TD	P9310M	P8310M	P7310M	P6310M	P5310M	P4310M	P3310M	 CN18, Direct (直连), 3m
GE-Ohmeda OxyTip sensor	P9110L	P8110L	P7110L	P6110L	P5110L	P4110L	P3110L	 CN17, GE Cable, 1m
S/5, i4 Modular Monitors with E-series modules; Corometrics 259 with TruSignal(Ohmeda tech) B20/30/40 Series	P9310K	P8310K	P7310K	P6310K	P5310K	P4310K	P3310K	 CN15, Direct (直连), 3m
<b>GE TruSignal</b>								
TruSignal sensor	P9110S	P8110S	P7110S	P6110S	P5110S	P4110S	P3110S	 CN64, GE Cable, 2m
	P9110T	P8110T	P7110T	P6110T	P5110T	P4110T	P3110T	 CN74, GE Cable, 3m
<b>Goldway 金科威</b>								
940 (Nellcor tech)	P9311	P8311	P7311	P6311	P5311	P4311	P3311	 CMS, Direct (直连), 3m
Nanjing SD, SENMIND 900F, KELIWEI, STRONG 7000C, HRD EM6/EM7, Snolife Smart wave15, CMS6500	P9311A	P8311A	P7311A	P6311A	P5311A	P4311A	P3311A	 CPS, Direct (直连), 3m
UT4000F (Digital)	P9311D	P8311D	P7311D	P6311D	P5311D	P4311D	P3311D	
(Nellcor tech)	P9311B	P8311B	P7311B	P6311B	P5311B	P4311B	P3311B	 CPT7, Direct (直连), 3m
UT4000B (Nellcor oximax)	P9311E	P8311E	P7311E	P6311E	P5311E	P4311E	P3311E	
<b>Guoteng 珠海国腾</b>								
Carewell, Nova-3D (BCI tech)	P9350	P8350	P7350	P6350	P5350	P4350	P3350	 CP6, Direct (直连), 3m
<b>hlmedical 海利康</b>								
TE-4000X	P91128	P81128	P71128	P61128	P51128	P41128	P31128	 CMS, Direct (直连), 1m

Reusable SpO2 Sensors/重复性血氧探头

Monitor 监护仪厂家	Adult clip 成人指夹	Adult soft 成人软指套	Pedia. Clip 儿童指夹	Pedia. Soft 儿童软指套	Neonate wrap 新生儿包裹	Multi-Y 多功能Y型	Adult Ear clip 成人耳夹	Connector/Adapter/Len. 接头/转接线/长度
<b>Hwatime 华腾</b>								
Hwatime UN-Medical galaxy	P9364A	P8364A	P7364A	P6364A	P5364A	P4364A	P3364A	 CP6R, Direct (直连), 3m
H6(BCI)	P93126	P83126	P73126	P63126	P53126	P43126	P33126	
<b>Infinium</b>								
OMNI 2, OMNI 3; STRONG(Digital)	P9365	P8365	P7365	P6365	P5365	P4365	P3365	 CP5R, Direct (直连), 3m
OMNI(BCI Tech)	P9365A	P8365A	P7365A	P6365A	P5365A	P4365A	P3365A	 CP6, Direct (直连), 3m
OMNI III(Digital)	P9365E	P8365E	P7365E	P6365E	P5365E	P4365E	P3365E	 CP5, Direct (直连), 3m
OMNI II	P9365B	P8365B	P7365B	P6365B	P5365B	P4365B	P3365B	 CN90, Direct (直连), 3m
OMNI III, Comen CS0/80 (Digital)	P9365C	P8365C	P7365C	P6365C	P5365C	P4365C	P3365C	 CN93, Direct (直连), 3m
<b>Innomed</b>								
Oxycard, Innocare-T0	P9369	P8369	P7369	P6369	P5369	P4369	P3369	 CP7R, Direct (直连), 3m
<b>Invivo 英威</b>								
4500, 4500 Plus, Millenia Omega 1445, Omni-Trak TVS, Sat Scout, 4500 Plus 2	P9312	P8312	P7312	P6312	P5312	P4312	P3312	 CM7, Direct (直连), 3m
New generation (Masim)	P9312A	P8312A	P7312A	P6312A	P5312A	P4312A	P3312A	 CM7, Direct (直连), 3m
<b>Ivy Biomedical</b>								
405P, 405T	P9319	P8319	P7319	P6319	P5319	P4319	P3319	 CN30, Direct (直连), 3m
Vital-Guard 450C	P9315	P8315	P7315	P6315	P5315	P4315	P3315	 CN22, Direct (直连), 3m
<b>Jeteem 井田医电</b>								
Bd9000	P9349	P8349	P7349	P6349	P5349	P4349	P3349	 CP6R, Direct (直连), 3m

Reusable SpO2 Sensors/重复性血氧探头

SpO2

Monitor 监护仪厂家	Adult clip 成人指夹	Adult soft 成人软指套	Pedia. Clip 儿童指夹	Pedia. Soft 儿童软指套	Neonate wrap 新生儿包裹	Multi-Y 多功能Y型	Adult Ear clip 成人耳夹	Connector/Adapter/Len. 接头/转接线/长度
<b>KERNEL</b>								
KERNEL	P9364B	P8364B	P7364B	P6364B	P5364B	P4364B	P3364B	 CP6, Direct (直连), 3m
<b>Konica Minolta 柯尼卡</b>								
Portable Pulsox-300, Pulsox-300i	P9177	P8177	P7177	P6177	P5177	P4177	P3177	 CN100, Direct (直连), 3m
<b>Konsung 江苏康尚</b>								
Aurora	P93133	P83133	P73133	P63133	P53133	P43133	P33133	 CN39, Direct (直连), 3m
<b>Kontron 康强</b>								
7840, 7845, Minimon 7137, Micromon, Supermon, Kolormon, Trakmon	P9313 0608010	P8313	P7313	P6313	P5313	P4313	P3313	 CN20, Direct (直连), 3m
7 pins 80 degree	P9313A	P8313A	P7313A	P6313A	P5313A	P4313A	P3313A	 CP7T, Direct (直连), 3m
<b>L&amp;T Medical</b>								
Stellar P, Star 50, Galaxy, Planet, Planet 50, Cosmic	P9376	P8376	P7376	P6376	P5376	P4376	P3376	 CN79, Direct (直连), 3m
	P9376A	P8376A	P7376A	P6376A	P5376A	P4376A	P3376A	
	P9376B	P8376B	P7376B	P6376B	P5376B	P4376B	P3376B	 CP6R, Direct (直连), 3m
<b>Masim 马西莫</b>								
GE PRO100Q, V300, NVIVO, TBI, Schiller, Datascope, ZOLL R Series, MDL Lifepak 20 (Masim tech)	P9315 1396/1969	P8315	P7315	P6315	P5315	P4315	P3315	 CN22, Direct (直连), 3m
(Masim tech)	P9115A	P8115A	P7115A	P6115A	P5115A	P4115A	P3115A	 CN23, See Masim cable, 1m
DB9PIN (Masim Oximax)	P9115F 1863/LNCS DC	P8115F	P7115F	P6115F	P5115F	P4115F	P3115F	
(Masim tech)	P9115B LNOP DC	P8115B	P7115B LNOP DC- IP	P6115B	P5115B	P4115B 1544	P3115B LNOP TC-I	 CN24, See Masim cable, 1m
(Masim Oximax)	P9115K	P8115K	P7115K	P6115K	P5115K	P4115K	P3115K	

Reusable SpO2 Sensors/重复性血氧探头

Reusable SpO2 Sensors/重复性血氧探头

SpO2

Monitor 监护仪厂家	Adult clip 成人夹夹	Adult soft 成人软指套	Pedia. Clip 儿童指夹	Pedia. Soft 儿童软指套	Neonate wrap 新生儿包裹	Multi-Y 多功能Y型	Adult Ear clip 成人耳夹	Connector/Adapter/Len. 接头/转接线/长度
<b>Masim 马西莫</b>								
(Masim Oximax)	P9115H	P8115H	P7115H	P6115H	P5115H	P4115H	P3115H	 CN86, See Masim cable, 1m
Radical-7, Rad-87, Rad-57, Pronto,兼容RED-8,ZOLL X Series (Masim Oximax)	P9315G	P8315G	P7315G	P6315G	P5315G	P4315G	P3315G	 CN80, Direct (直连), 3m
Corplus, Philips Efficia CM100, Zoll (Masim Oximax) (N)	P9315Q	P8315Q	P7315Q	P6315Q	P5315Q	P4315Q	P3315Q	 CN44, Direct (直连), 3m
RD rainbow, CM 100/120/150 Corplus, Philips Efficia CM100, Zoll (Masim Oximax) (N)	P9315Q-B	P8315Q-B	P7315Q-B	P6315Q-B	P5315Q-B	P4315Q-B	P3315Q-B	 CN44-B, Direct (直连), 3m
RD SET (Masim Oximax) (N)	P9115S RD SET DC-I	P8115S RD SET DB-I	P7115S RD SET DG-P	P6115S	P5115S	P4115S RD SET Y-4 Multi site Sensor	P3115S RD SET TC-I	 CN443, Direct (直连), 1m
Rad-G (N)	P9115AS	P8115AS	P7115AS	P6115AS	P5115AS	P4115AS	P3115AS	 CN443A, Direct (直连), 1m
Rad-97 (N)	P9315T	P8315T	P7315T	P6315T	P5315T	P4315T	P3315T	 CN457, Direct (直连), 3m
<b>Medlab</b>								
NANOX2, NANOXeco	P9303A	P8303A	P7303A	P6303A	P5303A	P4303A	P3303A	 CN81, Direct (直连), 3m
PEARL10, NANOX10, Nanox10C	P9319	P8319	P7319	P6319	P5319	P4319	P3319	 CN30, Direct (直连), 3m
<b>Medsorce</b>								
CAPNOX, VITRO, élan, POX 10, P-OX100L	P9303	P8303	P7303	P6303	P5303	P4303	P3303	 CP7, Direct (直连), 3m
BPM-200	P9319	P8319	P7319	P6319	P5319	P4319	P3319	 CN30, Direct (直连), 3m
<b>MEK</b>								
MP 100,110,400,1000	P9316	P8316	P7316	P6316	P5316	P4316	P3316	 CN30B, Direct (直连), 3m
	P9116	P8116	P7116	P6116	P5116	P4116	P3116	 CN30B, Direct (直连), 3m
Normal 8 pins	P9316A	P8316A	P7316A	P6316A	P5316A	P4316A	P3316A	 CN65, Direct (直连), 3m

Monitor 监护仪厂家	Adult clip 成人夹夹	Adult soft 成人软指套	Pedia. Clip 儿童指夹	Pedia. Soft 儿童软指套	Neonate wrap 新生儿包裹	Multi-Y 多功能Y型	Adult Ear clip 成人耳夹	Connector/Adapter/Len. 接头/转接线/长度
<b>MEK</b>								
Mp1000	P9316B	P8316B	P7316B	P6316B	P5316B	P4316B	P3316B	 CN65, Direct (直连), 3m
PAMO II	P9316C	P8316C	P7316C	P6316C	P5316C	P4316C	P3316C	 CN65, Direct (直连), 3m
<b>Mennen 美能</b>								
Horizon XL, Mercury 1200, MR 1300, Envoy	P9317	P8317	P7317	P6317	P5317	P4317	P3317	 CM10, Direct (直连), 3m
Normal 8 pins	P9317A	P8317A	P7317A	P6317A	P5317A	P4317A	P3317A	 CN25, Direct (直连), 3m
(Masim tech)	P9317B	P8317B	P7317B	P6317B	P5317B	P4317B	P3317B	 CN26, Direct (直连), 3m
(Oximax tech)	P9317D	P8317D	P7317D	P6317D	P5317D	P4317D	P3317D	
VitaLogjik (4000/4500/6000) (Masim tech)	P9317F	P8317F	P7317F	P6317F	P5317F	P4317F	P3317F	
<b>Mindray 迈瑞</b>								
12 Male pins	P9318	P8318	P7318	P6318	P5318	P4318	P3318	 CN94, Direct (直连), 3m
12 Male pins (Masim tech)	P9318A	P8318A	P7318A	P6318A	P5318A	P4318A	P3318A	
503,507,508,509 (Criticare tech)	P9318B	P8318B	P7318B	P6318B	P5318B	P4318B	P3318B	 CM5, Direct (直连), 3m
Normal 5 pins	P9318C	P8318C	P7318C	P6318C	P5318C	P4318C	P3318C	 CP5, Direct (直连), 3m
Normal 6 pins	P9318D	P8318D	P7318D	P6318D	P5318D	P4318D	P3318D	 CP6, Direct (直连), 3m
<b>Mindray 迈瑞</b>								
Mindray :PM7000, Pm8000, PM9000,MEC-1000,MEC-2000 Creative: UP-6000, Pc9000 Edan:M9,M99A,M8, M8A,M8B,M3,M3B	P9318E	P8318E	P7318E	P6318E	P5318E	P4318E	P3318E	 CP6R, Direct (直连), 3m
6 male pins, 40 degree (Masim tech)	P9318F	P8318F	P7318F	P6318F	P5318F	P4318F	P3318F	

Reusable SpO2 Sensors/重复性血氧探头

Monitor 监护仪厂家	Adult clip 成人指夹	Adult soft 成人软指套	Pedia. Clip 儿童指夹	Pedia. Soft 儿童软指套	Neonate 新生儿 新生儿包裹	Multi-Y 多功能Y型	Adult 成人耳夹	Connector/Adapter/Len. 接头/转接线/长度
<b>Mindray 迈瑞</b>								
6 male pins, 60 degree(Masim tech)	P9318G	P8318G	P7318G	P6318G	P5318G	P4318G	P3318G	 CP6S, Direct (直连), 3m
6 male pins, 80 degree(Nellcor Oximax)	P9318L	P8318L	P7318L	P6318L	P5318L	P4318L	P3318L	 CP6T, Direct (直连), 3m
PM-6800,BeneView T5 OEM No.0010-30-42737, VS600/VS900,iMEC/iMEC 12, Penlon SP M8,EDAN X12	P9318H	P8318H	P7318H	P6318H	P5318H	P4318H	P3318H	 CN27, Direct (直连), 3m
T5/T8(Masim tech)	P9318J	P8318J	P7318J	P6318J	P5318J	P4318J	P3318J	 CN28, Direct (直连), 3m
(Masim Oximax),DPM	P9318P	P8318P	P7318P	P6318P	P5318P	P4318P	P3318P	 CN9J, Direct (直连), 3m
T5/T8(Nellcor Oximax)	P9318K	P8318K	P7318K	P6318K	P5318K	P4318K	P3318K	 CN29, Direct (直连), 3m
<b>Mindray-Datascope</b>								
Accusat, Accutor 3/4/ SAT/PLUS Passport (Datascope Tech)	P9308	P8308	P7308	P6308	P5308	P4308	P3308	 CN07, Direct (直连), 3m
Accusat, Accutor 3/4/ SAT/PLUS Passport (Datascope Tech)	P9108A	P8108A	P7108A	P6108A	P5108A	P4108A	P3108A	 CN08, P0208, 1m
Passport II/IILT, XG, Accutor Plus, Expert(Nellcor) (Trio/Spectrum w/opt Nellcor)	P9319B	P8319B	P7319B	P6319B	P5319B	P4319B	P3319B	 CN22, Direct (直连), 3m
Passport II/IILT, XG, Accutor Plus, Trio, Spectrum (Masim)	P9315	P8315	P7315	P6315	P5315	P4315	P3315	 CN24, Direct (直连), 3m
Passport II/IILT, XG, Accutor Plus, Trio, Spectrum (Masim)	P9315B	P8315B	P7315B	P6315B	P5315B	P4315B	P3315B	 CN24, Direct (直连), 3m
<b>M&amp;B 北京麦邦</b>								
CD531, CD2000,MB526T	P9351	P8351	P7351	P6351	P5351	P4351	P3351	 CN7Q, Direct (直连), 3m
<b>Nellcor</b>								
NPB 40, NP820, NP875, NPB 180, NP8185, NPB-190, NPB-195, NPB 3900, NPB 4000	P9319	P8319	P7319	P6319	P5319	P4319	P3319	 CN30, Direct (直连), 3m
	P9119 512F-30- 28263	P8119 512E-30- 90391	P7119 512H-30- 79061	P6119	P5119 5188-30- 72107	P4119	P3119	 CN30, See cable, 1m

Reusable SpO2 Sensors/重复性血氧探头

SpO2

Monitor 监护仪厂家	Adult clip 成人指夹	Adult soft 成人软指套	Pedia. Clip 儿童指夹	Pedia. Soft 儿童软指套	Neonate 新生儿 新生儿包裹	Multi-Y 多功能Y型	Adult 成人耳夹	Connector/Adapter/Len. 接头/转接线/长度
<b>Nellcor</b>								
N5 50, N5 60, N5 95, N6 00, N7 51, N5 500 DASH 2 500, SOLAR 8 000, 8 000M, N6 5 Except CARESCAPE V1 00	P9319A	P8319A	P7319A	P6319A	P5319A	P4319A	P3319A	 CN31, Direct (直连), 3m
	P9119A DS-100A	P8119A	P7119A	P6119A	P5119A	P4119A D-YS	P3119A	 CN31, See Nellcor, 1m
NPB3000-3100, Symphony, NPB290,295,395	P9319B	P8319B	P7319B	P6319B	P5319B	P4319B	P3319B	 CN22, Direct (直连), 3m
DOC-10 direct sensor Mediana YM-2002 Series	P9319C	P8319C	P7319C	P6319C	P5319C	P4319C	P3319C	 CN46, Direct (直连), 3m
<b>Neusoft 东软</b>								
Nova-3D	P9343	P8343	P7343	P6343	P5343	P4343	P3343	 CP7, Direct (直连), 3m
General(Digital) Newtech NT3A	P9343A	P8343A	P7343A	P6343A	P5343A	P4343A	P3343A	 CN39 Direct (直连), 3m
(Nellcor Oximax)	P9343B	P8343B	P7343B	P6343B	P5343B	P4343B	P3343B	 CP8, Direct (直连), 3m
(Digital)	P9343C	P8343C	P7343C	P6343C	P5343C	P4343C	P3343C	 CP8, Direct (直连), 3m
<b>Nihon Kohden 光电</b>								
AL-800PA Module; BSM-1101, 2101/2, 3101,7100, 7200; OLV-1100A, Lifescope	P9321 TL-101S	P8321	P7321	P6321	P5321	P4321	P3321	 CN32, Direct (直连), 3m
BSM-2301K, BSM-2303K, BSM-2351, BSM-2353, BSM-4103, BSM-4113, VISMO, PVM-2701	P9121A TL-201T	P8121A	P7121A	P6121A	P5121A TL-220T	P4121A	P3121A	 CN33, P0221B, 1m
	P9321B	P8321B	P7321B	P6321B	P5321B	P4321B	P3321B	 CN34, Direct (直连), 3m
Lifescope OPV 1500 1500k	P9321C	P8321C	P7321C	P6321C	P5321C	P4321C	P3321C	 CN84, Direct (直连), 3m
BSM Series (Nellcor Oximax)	P9321D	P8321D	P7321D	P6321D	P5321D	P4321D	P3321D	 CN85, Direct (直连), 3m
BSM Series (Masim)	P9321E	P8321E	P7321E	P6321E	P5321E	P4321E	P3321E	 CN85, Direct (直连), 3m



Reusable SpO2 Sensors/重复性血氧探头

Monitor 监护仪厂家	Adult dip 成人探头	Adult soft 成人软探头	Pedia. Clip 儿童夹	Pedia. Soft 儿童软探头	Neonate wrap 新生儿包裹	Multi-Y 多功能Y型	Adult Ear clip 成人耳夹	Connector/Adapter/Len. 接头/转接线/长度
<b>Nihon Kohden 光电</b>								
SVM-7601	P9321F	P8321F	P7321F	P6321F	P5321F	P4321F	P3321F	CP6R, Direct (直连), 3m
<b>OXIMAX</b>	P9321G	P8321G	P7321G	P6321G	P5321G	P4321G	P3321G	CN22, Direct (直连), 3m
<b>Nonin 燕牌</b>								
3150	P9022C	P8022C	P7022C	P6022C	P5022C	P4022C	P3022C	CN54, Direct (直连), 0.35m
8604	P9322	P8322	P7322	P6322	P5322	P4322	P3322	CM6, Direct (直连), 3m
8500, PamSAT 2500, Avant 9600, Avant 9700, Avant 2120, 9840	P9122B 8000AA-1+85.00M	P8122B	P7122B	P6122B	P5122B	P4122B	P3122B	CN83, Direct (直连), 1m
<b>Northern Meditech 诺然美善</b>								
Northern Meditech: Gemini, Taurus, Virgo, Aquarius Progetti (Italy): pg S10 aXcent (Germany): Cetus X12	P9392	P8392	P7392	P6392	P5392	P4392	P3392	CN98, Direct (直连), 3m
<b>Novamatrix 诺医</b>								
505, 510, 511, 515B, 515C, 520, 7000, 7100, Oxypleth	P9323 8776-00	P8323	P7323	P6323	P5323	P4323	P3323	CN36, Direct (直连), 3m
Dixtal (Brazil)	P9323C	P8323C	P7323C	P6323C	P5323C	P4323C	P3323C	CN36B, Direct (直连), 3m
512/513/710/715, MARS	P9123A	P8123A	P7123A	P6123A	P5123A	P4123A	P3123A	CN37, Direct (直连), 1m
512/513/710/715, MARS	P9123B	P8123B	P7123B	P6123B	P5123B	P4123B	P3123B	CN04, P0223, 1m
<b>Osen 深圳奥生</b>								
Osen, Carewell CPM-9000T	P9359	P8359	P7359	P6359	P5359	P4359	P3359	CP6R, Direct (直连), 3m
Osen	P9359A	P8359A	P7359A	P6359A	P5359A	P4359A	P3359A	CP6R, Direct (直连), 3m
9000D	P9359B	P8359B	P7359B	P6359B	P5359B	P4359B	P3359B	CP6, Direct (直连), 3m

Reusable SpO2 Sensors/重复性血氧探头

SpO2

Monitor 监护仪厂家	Adult dip 成人探头	Adult soft 成人软探头	Pedia. Clip 儿童夹	Pedia. Soft 儿童软探头	Neonate wrap 新生儿包裹	Multi-Y 多功能Y型	Adult Ear clip 成人耳夹	Connector/Adapter/Len. 接头/转接线/长度
<b>Pace Tech Generra</b>								
Vitalmax 800 Series	P9311C	P8311C	P7311C	P6311C	P5311C	P4311C	P3311C	CN19, Direct (直连), 3m
<b>Palco 派克</b>								
All	P9324	P8324	P7324	P6324	P5324	P4324	P3324	CN38, Direct (直连), 3m
Palco	P9324A	P8324A	P7324A	P6324A	P5324A	P4324A	P3324A	CN67, Direct (直连), 3m
<b>Philips-HP- Agilent 飞利浦</b>								
78352C, 78354C, 78834C, M1020A, M1025B, M1025A, M1350B, M1722A/B, M1732A/B, M2475B, ViridaA1, A3	P9325 M1900B+ M11967	P8325 M1900B+ M1191T	P7325	P6325 M1900B+ M1192T	P5325 M1900B+ M1193T	P4325	P3325	CN39, Direct (直连), 3m
MP20, MP30, MP40/50/60/70/90, VM4, VM6, VM8, M3046A	P9325A M1196A	P8325A M1191A/B M1191A/L	P7325A M1196A	P6325A M1192A	P5325A M1193A	P4325A	P3325A M1194A	CN40, Direct (直连), 3m
	P9125A	P8125A	P7125A	P6125A	P5125A	P4125A	P3125A	CN40, Direct (直连), 1m
Philips-Masim	P9325B	P8325B	P7325B	P6325B	P5325B	P4325B	P3325B	CN41, Direct (直连), 3m
	P9325L	P8325L	P7325L	P6325L	P5325L	P4325L	P3325L	CN41, Direct (直连), 3m
<b>PETAS</b>								
KMA 800	P9379B	P8379B	P7379B	P6379B	P5379B	P4379B	P3379B	CP6R, Direct (直连), 3m
<b>Physio-Control</b>								
Lifepak 12	P9319	P8319	P7319	P6319	P5319	P4319	P3319	CN30, Direct (直连), 3m
Lifepak 12, 20 (Masim tech)	P9315	P8315	P7315	P6315	P5315	P4315	P3315	CN22, Direct (直连), 3m
<b>Perlong 普朗</b>								
EDWNING EMP00000 (澳威) UN-Medical galaxy (威达) Medcon DK-60005 (曼德隆) MILLION ML-1200 (美辰) IM 12E (亿康) (韩国) KP-9000 (韩国) Sunray FO-M99A (三能) P51-3000, S2000 P51-12, S2000 PRO-META, Perlong P51-3000, Zonare, Uric, Unicare, Beijing, Henry	P9364	P8364	P7364	P6364	P5364	P4364	P3364	CP6, Direct (直连), 3m
<b>Primedic 普美康</b>								
Primedic	P9375	P8375	P7375	P6375	P5375	P4375	P3375	CP6R, Direct (直连), 3m

Reusable SpO2 Sensors/重复性血氧探头

Monitor 监护仪厂家	Adult clip 成人指夹	Adult soft 成人软指套	Pedia. Clip 儿童指夹	Pedia. Soft 儿童软指套	Neonate wrap 新生儿包裹	Multi-Y 多功能Y型	Adult Ear clip 成人耳夹	Connector/Adapter/Len. 接头/转接线/长度
<b>Ruibo (Hunan) 湖南瑞博</b>								
16#12.5mm (BCI tech)	P9367	P8367	P7367	P6367	P5367	P4367	P3367	 CN68, Direct (直连), 3m
12#9mm, PM-9000A+ (BCI tech)	P9367A	P8367A	P7367A	P6367A	P5367A	P4367A	P3367A	 CN73, Direct (直连), 3m
12#9mm, PM-9000A+ (Digital)	P9367B	P8367B	P7367B	P6367B	P5367B	P4367B	P3367B	 CN40, Direct (直连), 3m
<b>Schiller 摩轴</b>								
RGB Omicrom FT, Omicron PREMIUM	P9303	P8303	P7303	P6303	P5303	P4303	P3303	 CP7, Direct (直连), 3m
Schiller Truscope (Nelcor oximax)	P9326C	P8326C	P7326C	P6326C	P5326C	P4326C	P3326C	
CM8,Schiller,Argus, TM-7(Nelcor tech)	P9326	P8326	P7326	P6326	P5326	P4326	P3326	 CN42, Direct (直连), 3m
TM7,Argus LCM (Masim)	P9326A	P8326A	P7326A	P6326A	P5326A	P4326A	P3326A	 CN43, Direct (直连), 3m
TM 910	P9326B	P8326B	P7326B	P6326B	P5326B	P4326B	P3326B	
<b>Spacelabs 太空</b>								
90419, 90501,90651A-0390351-0/6, 90465, 90466, 90467,90651A-08, PC Express, 90489, 1m77	P9327 015-0130-00+700-0030-00	P8327	P7327	P6327	P5327	P4327	P3327	 CN44, Direct (直连), 3m
90496 Ultraview  (Maximo non-oximax)	P9327A 015-0130-00+700-0030-00	P8327A	P7327A	P6327A	P5327A	P4327A	P3327A	 CN45, Direct (直连), 3m
Ultraview SL (Maximo Oximax)	P9327E	P8327E	P7327E	P6327E	P5327E	P4327E	P3327E	 CN69, Direct (直连), 3m
Mcare300 (Oximax tech)	P9327D	P8327D	P7327D	P6327D	P5327D	P4327D	P3327D	 CN72, Direct (直连), 3m
<b>Solaris 索莱瑞</b>								
(Digital)	P9320	P8320	P7320	P6320	P5320	P4320	P3320	 CP5, Direct (直连), 3m

Reusable SpO2 Sensors/重复性血氧探头

SpO2

Monitor 监护仪厂家	Adult dip 成人指夹	Adult soft 成人软指套	Pedia. Clip 儿童指夹	Pedia. Soft 儿童软指套	Neonate wrap 新生儿包裹	Multi-Y 多功能Y型	Adult Ear clip 成人耳夹	Connector/Adapter/Len. 接头/转接线/长度
<b>Solaris 索莱瑞</b>								
Solaris 索莱瑞	P9320F	P8320F	P7320F	P6320F	P5320F	P4320F	P3320F	 CP5, Direct (直连), 3m
General G2A, G3A, G3C, G3H, G3F (Digital)	P9320A	P8320A	P7320A	P6320A	P5320A	P4320A	P3320A	 CP6R, Direct (直连), 3m
(Digital)	P9120B	P8120B	P7120B	P6120B	P5120B	P4120B	P3120B	 CN30, Direct (直连), 1m
Spurt (Digital)	P9320C	P8320C	P7320C	P6320C	P5320C	P4320C	P3320C	 CP7, Direct (直连), 3m
Creative UP-7000, PC-5000A (Digital)	P9320D	P8320D	P7320D	P6320D	P5320D	P4320D	P3320D	 CP5T, Direct (直连), 3m
(Digital)	P9320E	P8320E	P7320E	P6320E	P5320E	P4320E	P3320E	 CN40, Direct (直连), 3m
<b>Senmind 深迈</b>								
KP-9000	P9348	P8348	P7348	P6348	P5348	P4348	P3348	 CN40, Direct (直连), 3m
(Digital)	P9348A	P8348A	P7348A	P6348A	P5348A	P4348A	P3348A	
	P9348B	P8348B	P7348B	P6348B	P5348B	P4348B	P3348B	 CP6, Direct (直连), 3m
<b>Tianrong Medical 徐州天荣</b>								
TR-900E	P9357	P8357	P7357	P6357	P5357	P4357	P3357	 CP5, Direct (直连), 3m
Tianrong Medical	P9357A	P8357A	P7357A	P6357A	P5357A	P4357A	P3357A	 CP6R, Direct (直连), 3m
<b>UTECH 重庆如泰</b>								
UTECH	P93127	P83127	P73127	P63127	P53127	P43127	P33127	 CP6, Direct (直连), 3m
<b>Welch Allyn Protocol 伟伦</b>								
Propaq 206-EL, 204-EL, 202-EL, 106-EL	P9319	P8319	P7319	P6319	P5319	P4319	P3319	 CN30, Direct (直连), 3m

Reusable SpO2 Sensors/重复性血氧探头

Monitor 监护仪厂家	Adult clip 成人指夹	Adult soft 成人软指套	Pedia. Clip 儿童指夹	Pedia. Soft 儿童软指套	Neonate wrap 新生儿包裹	Multi-V 多功能V型	Adult Ear clip 成人耳夹	Connector/Adapter /Len. 接头/转接线/长度
<b>wisonic 华声</b>								
wisonic	P91104	P81104	P71104	P61104	P51104	P41104	P31104	 CN92, Direct (直连), 1m
(Digital)	P93104	P83104	P73104	P63104	P53104	P43104	P33104	 CN92, Direct (直连), 3m
<b>Worldlife</b>								
Worldlife	P93134	P83134	P73134	P63134	P53134	P43134	P33134	 CN495, Direct (直连), 3m
<b>Yonker 永康</b>								
Waltvick Zafiro	P93124	P83124	P73124	P63124	P53124	P43124	P33124	 CN438, Direct (直连), 3m
<b>Zoncare 中旗</b>								
(BCI)	P9370	P8370	P7370	P6370	P5370	P4370	P3370	 CM6-OB, Direct (直连), 3m
(BCI)	P9370A	P8370A	P7370A	P6370A	P5370A	P4370A	P3370A	 CN92, Direct (直连), 3m
Zoncare PM-7000M (Digital)	P9370B	P8370B	P7370B	P6370B	P5370B	P4370B	P3370B	 CN92, Direct (直连), 3m
<b>Zondan 和心重典</b>								
Zd120	P9345	P8345	P7345	P6345	P5345	P4345	P3345	 CP6, Direct (直连), 3m
iMD12	P9345B	P8345B	P7345B	P6345B	P5345B	P4345B	P3345B	 CN89, Direct (直连), 3m